

# Scaffold Analysis of Isobaric Tag Datasets with a Pooled Reference

Open Scaffold, choose Edit/Preferences/Memory, then set Memory to 50-80% of RAM, with a recommended setting of  $\geq 4000$  MB. Close, re-open Scaffold and then open the Scaffold file. Choose Protein Threshold: 99%, Min. # Peptides: 2, Peptide Threshold: 95%



Click on Q+, enter these settings after clicking Quant/Update Experimental Design: Analysis Type: Intensity Based Experiment: Between Subjects with (Common/Pooled Reference)



## “Edit Sample Names and Categories”

Select the number of categories or conditions in the experiment *including* the pooled reference. Double-click on each cell to rename each quant sample under the “Name” column, which corresponds to its TMT label



## “Organize Quant Samples”

Move Reference, Control, & Experimental Samples into proper categories For experiments with multiple mastermixes be sure to click “Add reference alignment group” for each new mastermix and corresponding samples



## “Approve Settings”

Click “Edit Settings”, on the Normalization Tab choose Calculation Type: Median, Blocking Level: Unique Peptides, Reference Type: Average Protein Reference, Spectrum Quality Filter: Reference value required, Normalization Between Samples: On. On the Minimum Dynamic Range tab chose 5%; on the Other Settings tab deselect “Use Non-Exclusive Peptides” to *exclude* peptides that below to multiple proteins. Click “Apply” and “Finish”



## “Quant/Quantitative Testing”

Select two groups of samples to compare (e.g., treated vs control), use the “Remove” button to remove other samples. Choose Mann-Whitney Test with Benjamini-Hochberg multiple test correction with a significance level of  $p < 0.05$ , then choose the appropriate reference Category (e.g., control) and click “Apply”.



## Interpreting Results

Set Display Options to Fold Change Ratio, Min FC: 1.5 Fold; Click on Mann-Whitney Test column to rank order proteins by p-value. In addition to the Mann-Whitney test, Scaffold also carries out a Brown-Forsythe population variance test. Orange highlighting indicates proteins whose p-values have significantly *unequal* variances based on this test. The Scaffold tutorial cautions that these p-values “may be based on a small number of measurements or unreliable data” and that they should be examined in more depth.