

***Analysis of Protein Complexes
by Size Exclusion
Chromatography Coupled with
Light Scattering***


Size Exclusion Chromatography (SEC) Coupled with Light Scattering (LS)

- *Derivation of Molecular Weight from LS experiment*
- *Experimental Set Up for SEC/LS “in-line”*
- *Evaluation of the SEC/LS System*
 - Results for Standard Proteins*
 - Sample Requirements*
- *Applications of SEC/LS to study protein complexes*
- *Conclusions*

Size Exclusion Chromatography (SEC) Coupled with Light Scattering (LS)

- ***Derivation of Molecular Weight from LS experiment***
- *Experimental Set Up for SEC/LS “in-line”*
- *Evaluation of the SEC/LS System*
 - Results for Standard Proteins*
 - Sample Requirements*
- *Applications of SEC/LS to study protein complexes*
- *Conclusions*

Light Scattering Experiments



Static Light Scattering Experiments

Debye-Zimm formalism for $R(\Theta)$, the excess intensity of scattered light at an angle Θ

$$\frac{K^*c}{R(\theta)} = \frac{1}{M_w P(\theta)} + 2A_2c$$

c is the sample concentration (g/ml)

M_w is the weight-average molecular weight (molar mass)

A_2 is the second virial coefficient (ml-mol/g²)

K^* is an optical parameter equal to $4\pi^2 n^2 (dn/dc)^2 / (\lambda_0^4 N_A)$

n is the solvent refractive index and dn/dc is the refractive index increment

N_A is Avogadro's number

λ_0 is the wavelength of the scattered light in vacuum (cm)

$P(\theta)$ is the form factor (describes angular dependence of scattered light)

Static Light Scattering Experiments

at low concentrations $c < 0.1$ mg/mL

$$2A_2cM_w \ll 1$$

thus, the second virial coefficient term ($2A_2c$) can be neglected

$$\frac{K^*c}{R(\theta)} = \frac{1}{M_w P(\theta)} + 2A_2c \quad \longrightarrow \quad \frac{K^*c}{R(\theta)} = \frac{1}{M_w P(\theta)}$$

expansion of $P(\Theta)$ to the first order gives

$$1/P(\Theta) = 1 + (16\pi^2/3\lambda^2) \langle r_g^2 \rangle \sin^2(\Theta/2) + \dots$$

Static Light Scattering Experiments

$$\frac{K^*c}{R(\theta)} = \frac{1}{M_w} \left(1 + \left(\frac{16\pi^2}{3\lambda^2} \right) \langle r_g^2 \rangle \sin^2\left(\frac{\theta}{2}\right) \right)$$


Using a multi angle instrument
construct a plot of

$$\frac{K^*c}{R(\theta)} \quad \text{against} \quad \sin^2\left(\frac{\theta}{2}\right)$$

From intercept \longrightarrow ***Derived MW***

Zimm Plot Ovalbumin (43 kDa)

$$\frac{K^*c}{R(\theta)} = \frac{1}{M_w} (1 + f(\sin^2(\frac{\theta}{2})))$$



At low concentrations

$$\frac{K * c}{R(\theta)} \quad \text{against} \quad \sin^2\left(\frac{\theta}{2}\right)$$

From intercept \longrightarrow ***Derived MW***


weight-average MW



fractionate samples


Size Exclusion Chromatography (SEC) Coupled with Light Scattering (LS)

- *Derivation of Molecular Weight from LS experiment*
- ***Experimental Set Up for SEC/LS “in-line”***
- *Evaluation of the SEC/LS System*
 - Results for Standard Proteins*
 - Sample Requirements*
- *Applications of SEC/LS to study protein complexes*
- *Conclusions*



Zimm Plot Ovalbumin (43 kDa)

$$\frac{K^*c}{R(\theta)} = \frac{1}{M_w} (1 + f(\sin^2(\frac{\theta}{2})))$$




Volume : 16.300 mL

Conc. : (0.173 ± 0.000) mg/mL

Mw $(42.79 \pm 0.03) \times 10^3$ g/mol


Radius : 0.0 ± 0.0 nm



Molar Mass Distribution Plot

Ovalbumin 43 kDa


Molar Mass vs. Volume



Molar Mass Distribution Plot

BSA 66 kDa

Molar Mass vs. Volume



Size Exclusion Chromatography (SEC) Coupled with Light Scattering (LS)

- *Derivation of Molecular Weight from LS experiment*
- *Experimental Set Up for SEC/LS “in-line”*
- ***Evaluation of the SEC/LS System***
 - Results for Standard Proteins*
 - Sample Requirements*
- *Applications of SEC/LS to study protein complexes*
- *Conclusions*

**Molecular Weights Determined from "in line" analyses;
static LS with SEC in line; 16 protein standards, MW 6.5 to 475 kDa**

Protein	Oligomeric state	# Runs	Pred. MW (kDa) ^a	Average MW \pm St. Dev. (kDa)	Average error (%)
Aprotinin	monomer	2	6.5	6.8 \pm 0.5	4.6
Cytochrome C	monomer	5	12.3	12.01 \pm 0.57	2.4
α -Lactalbumin	monomer	2	14.2	14.32 \pm 0.01	0.9
Myoglobin	monomer	3	17.0	14.19 \pm 0.91	16
β Lactglobulin	monomer	2	18.3	20.06 \pm 0.33	9.7
Tripsin inhibitor	monomer	1	20.0	20.50	2.3
Carbonic anhydrase	monomer	4	29.0	29.22 \pm 0.20	0.8
Ovalbumin	monomer	10	42.8	42.52 \pm 0.68	1.4
BSA (monomer)	monomer	5	66.4	66.41 \pm 1.00	1.2
Transferrin	monomer	2	75.2	76.92 \pm 0.98	2.3
Enolase (yeast)	dimer	3	93.3	80.74 \pm 1.18	13
Enolase (rabbit)	dimer	4	93.7	86.44 \pm 1.90	7.8
BSA (dimer)	dimer	5	132.9	137.10 \pm 3.93	3.2
Alc. dehydrogenase	tetramer	4	147.4	144.02 \pm 0.86	2.4
Aldolase (rabbit)	tetramer	2	156.8	153.7 \pm 1.91	1.1
Apo-ferritin	24 ^x monomer	2	475.9	470.3 \pm 2.62	1.2
Median error:					2.3

Buffer: 20 mM HEPES, 150 mM KCl, 1 mM EDTA, pH=8.0; column: Superdex 200 or Superdex 75

Correlation between the amount of protein analyzed and the accuracy of MW determination

Protein	Amount loaded (µg)	# Runs	Pred. MW (kDa)	Avrg. MW (kDa)	SD (kDa)	Avrg. error (%)	Range of accuracy (%)
Ovalbumin	150	4	42.8	42.4	0.3	0.9	0.2 to 1.6
	100	7	42.8	42.3	0.8	1.2	0.2 to 2.4
	45-50	4	42.8	41.6	1	2.8	0.5 to 5.8
	6-10	5	42.8	42.9	2	0.2	1.4 to 4.5
Transferrin	100	3	75.2	76.5	1	1.7	0.7 to 3.2
	8	5	75.2	76.3	2	1.5	0.3 to 5.2

column: TSK GEL G3000_{SWXL} [TosoHaas], buffer: 20 mM phosphate, 150 mM NaCl, pH=7.5

Sample requirements for proteins.

Column	Optimal amount of protein			Total volume of the eluting peak
	for expected MW >40 kDa	for expected MW 10 - 40 kDa	for expected MW <10 kDa	
Superose 6 (Pharmacia)	100 µg	N/A	N/A	~ 2mL
Superdex 200 (Pharmacia)	100 µg	200 - 300 µg	N/A	~ 2mL
Superdex 75 (Pharmacia)	50 µg	100 - 200 µg	400 µg	~ 1mL

Size Exclusion Chromatography (SEC) Coupled with Light Scattering (LS)

- *Derivation of Molecular Weight from LS experiment*
- *Experimental Set Up for SEC/LS “in-line”*
- *Evaluation of the SEC/LS System*
 - Results for Standard Proteins*
 - Sample Requirements*
- ***Applications of SEC/LS to study protein complexes***
- *Conclusions*

Applications of SEC/LS to study protein complexes


- *Determination of the oligomeric state of mutant vs. wild type protein*
- *Mixtures of non-interacting proteins*
- *Mixtures of interacting protein- detection of ligand driven protein complexes*
- *Determination of oligomeric state of glycosylated proteins*
- *Determination of oligomeric state of membrane proteins solubilized in detergents*

Determination of the oligomeric state of mutant vs. wild type protein

Example:

protein 12 kDa (WT protein exists as a trimer)

Three mutants and WT protein were analyzed.



Mixtures of non-interacting proteins

Example:

BSA monomer - 66 kDa protein


*Yeast Enolase - 93 kDa dimer
(2x46kDa)*

Analysis of co-eluting protein mixture

BSA 66kDa

BSA+ENY mixture

ENY dimer 93 kDa



Mixtures of interacting protein- detection of ligand driven protein oligomerization


Example:

*protein 27 kDa (protein exists as a
mixture of monomer and dimer)*


ligand 7 kDa

*Ligand binding shifts the protein into
dimeric form*


Analysis of interacting proteins




Molar Mass vs. Volume




Molar Mass vs. Volume



Molar Mass vs. Volume



Molar Mass vs. Volume



Determination of the oligomeric state of glycosylated protein

Data Analysis:

Use “three detector method”

Use ASTRA

(knowing the amount of sugars bound)

use weight-average dn/dc value

Three Detector Method

Jie Wen, Tsutomu Arakawa and John S. Philo

Anal Biochem 1996 Sep 5;240(2):155-66


Yutaro Hayashi, Hideo Matsui and Toshio Takagi

Methods Enzymol 1989;172:514-28

$$M_p = \frac{k^*(LS)(UV)}{\varepsilon(RI)^2}$$

$$MW_p = 91.39 \times [(LS) \cdot (UV) / (A \cdot (RI^2))]$$

Protein	MW (kDa)
Ova	43
BSA(1)	66
BSA(2)	132
Ald	156
Apo-Fer	475




Determination of the oligomeric state of glycosylated protein

Example:

<i>protein</i>	<i>58 kDa</i>	<i>extracellular ANP-binding domain (ECD) of cell-surface receptor</i>
		<i>16% of mass is sugar</i>
		<i>$dn/dc_t = 0.179 \text{ g/mL}$</i>
	<i>48 kDa</i>	<i>polypeptide portion</i>
<i>ligand</i>	<i>2.7 kDa</i>	<i>atrial natriuretic peptide (ANP)</i>


Injected sample complex (ECD : ANP) 2:1



$ECD_{dimer} = 2 \times 58 = 116 \text{ kDa}$ (96 kDa)

$ANP = 2.7 \text{ kDa}$

Injected sample complex (ECD : ANP) 2:1



$$ECD_{dimer} = 2 \times 58 = 116 \text{ kDa (96 kDa)}$$

$$ANP = 2.7 \text{ kDa}$$

ECD-ANP complex; ~3 μ g


MW_{glycoprotein} = 103 \pm 10 kDa

MW_{polypeptide} = 96 \pm 7 kDa

ECD; ~2 μ g

MW_{glycoprotein} = 54 \pm 6 kDa

MW_{polypeptide} = 44 \pm 5 kDa



ECD_{dimer} = 2 x 58 = 116 kDa (96 kDa)

ANP = 2.7 kDa

Determination of the oligomeric state of detergent solubilized membrane protein

Data Analysis:

Use “three detector method”

Use ASTRA

use “corrected” dn/dc value as described by Habayashi (scaled RI signal such that it represents contribution only from polypeptide)

Three Detector Method

Yutaro Hayashi, Hideo Matsui and Toshio Takagi

Methods Enzymol 1989;172:514-28

$$M_p = \frac{k^*(LS)(UV)}{\varepsilon(RI)^2}$$

allows determination of MW polypeptide
(oligomeric state of the protein)

allows determination of mass of detergent/lipids bound to polypeptide

$$\left(\frac{dn}{dc}\right)_{corr} = \left(\frac{dn}{dc}\right)_{pp} + \delta \left(\frac{dn}{dc}\right)_{d+l} = K \varepsilon \frac{(RI)}{(UV)}$$

δ is mass of detergent and/or lipids per 1 gram of polypeptide

Assumption : detergent does not produce any signal in UV

Determination of the oligomeric state of detergent solubilized protein

Example:

protein *47 kDa* *well characterized porin*

detergent


dodecyl maltoside (C12M) MW = 511 g/mol

0.5g/L i.e. 0.05%

CMC = 0.008% micelle size 50-70 kDa

porin monomer = 47 kDa

MW = 149 ± 3 kDa trimer




— *LS @ 90 degree*


— *RI*

— *UV @ 280 nm*

Peak ID - PRNC__DC



C12M micelle size = 72 kDa



Determination of the oligomeric state of detergent solubilized protein

Example:

protein *33 kDa*

Detergent

dodecyl maltoside (C12M) MW = 511 g/mol


n-Dodecyl- β -D-Maltoside


0.5g/L i.e. 0.05%

CMC = 0.008% micelle size 50-70 kDa


monomer = 33 kDa

MW = 225 ± 13 kDa heptamer





δ is mass of detergent and/or lipids per 1 gram of polypeptide



$MW_{complex} = 297 \text{ kDa}$

$MW_{polypeptide} = 149 \text{ kDa}$


$\delta = 1.02$

$MW_{complex} = 283 \text{ kDa}$

$MW_{polypeptide} = 225 \text{ kDa}$

$\delta = 0.26$

Molar Mass vs. Volume



Conclusions

SEC coupled with Static LS/RI/UV

- fast and accurate determination of molecular weight (MW) of macromolecules in solution
- single SEC/LS measurement should be sufficient to estimate a MW with a precession of $\pm 5\%$
- SEC/LS suitable for detection and characterization of non-interacting and interacting systems
- SEC/LS/UV/RI analysis can determine oligomeric state of detergent solubilized membrane proteins

Ken Williams

*Director of HHMI Biopolymer & W.M. Keck Biotechnology Resource
Laboratory*

NIH

Users of SEC/LS Service

Wyatt Technology